Article #00901S (US)
Produit #20901 (Canada)


The SNS™

Supplemental Nursing System For Breastfeeding Assistance
Sistema de lactancia complementario para dar el pecho
Système d'allaitement supplémentaire

Instructions for Use

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Mode d'emploi

■■■medela

TO THE NEW MOTHER

Dear Mother,

Congratulations! We at Medela wish many years of health and happiness for you and your baby. We are happy that you have chosen to breastfeed your baby and are pleased to be of help in providing the best nutrition for your baby. The Medela SNS™ breastfeeding assistance kit can help in nursing your adopted baby or in overcoming temporary breastfeeding difficulties, and will allow you and your baby to enjoy that special feeling of love, warmth, and trust that makes nursing so special. Please read the following instructions for use carefully.

Caution

You are using your SNS to help assist you in breastfeeding your baby. Baby should easily be making weight gains each week. If for some reason, your baby is not gaining weight, or if baby loses weight, consult with your health-care professional immediately. Consult with your health-care professional throughout the time you use your SNS to determine how often your baby should be checked and perhaps weighed. Please make sure you read this manual well, for it contains important information you need to know.

Your Doctor

Your physician and other health care providers should be made aware that you are using the SNS. Consult with them before you use the SNS as this will enable them to provide you with the best possible care. Follow your physician's recommendation as to the type of supplement to use, how much to use, what size tubing to use and any other suggestions they may have.

Table of Contents

ENGLISH

To the New Mother	1
Your SNS™	2
Additional Ordering Information	3
Getting Started	4-7
Cleaning and Sterilization	7
Specific Problems	12
Decreasing Supplements	16
Breastpumps	17

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
YOUR SNS

The SNS is designed to stimulate the mother's milk supply by providing the baby with the motivation needed to suck well.

Your SNS should include the following items:

1. Neck Cord
2. Feeding bottle
3. Three sizes of tubing
4. Valve holder
5. Attachment ring
6. Cover
7. Two rolls of hypoallergenic papertape
8. Instructions for use

ADDITIONAL ORDERING INFORMATION

To order parts for your SNS contact your local authorized Medela Rental or Retail Station.

<u>Part #</u>	<u>Item</u>
00901S(US)	Complete SNS
20901(Canada)	Complete SNS
8000794	Container with lid
8000041	Valve (white) with medium tubing
8000042	Valve (red) with small tubing
8000044	Valve (clear) with large tubing
8100084	Valve holder
8100085	Attachment ring
6000155	Neck Cord
8100086	Cover
8990007	Hypoallergenic papertape 1/2" x 10 meter
1907208	Instructions


Parts may also be ordered by calling 1-800-435-8316 and ordering with Discover Card, Mastercard, Visa, or American Express.

*Call customer service for current pricing
and shipping charges.*


GETTING STARTED

Before using the SNS, it may help if you get acquainted with the device. Practice assembling and disassembling your SNS as this will give you familiarity with the system.


1 Connect proper sized tubing (see Tubing Size). Generally, the middle size will be most appropriate. The smaller size is used for a slower flow or a thinner fluid. Consult with your health care professional for your proper size.

2 Push valve with attached tubings over valve holder. **Caution: Do not pull on tubing.** Tubings are not very strong and could break.


3 Slide tubing through attachment ring and push valve holder with the valve into the ring until it snaps.

4 Fill bottle with supplement recommended by physician. There is normally no need to warm the supplement. The tubing is very thin. Your body will heat the supplement as it passes through the tubing because the tubing is close to your skin. When taken from the

refrigerator or when supplement was stored in cold weather, warm the supplement to approximately body temperature.

- 5 Screw attachment ring onto feeding bottle, clamp both tubings into attachment ring to avoid milk flow as you place system around your neck. To avoid stretching make sure that the tubing has a loop.
- 6 Hang the SNS around your neck.

7 Place tubing on the top of the nipple making sure it extends about 1/4 inch beyond nipple and secure vertically with a piece of tape down to the base of the areola. Another piece of tape may be placed horizontally at base of areola. Put tape on breast and areola but not on the nipple itself. Most babies are not bothered by the tape.

8 Position baby at breast and release just the tubing on the breast you are about to nurse. Place all your fingers underneath the breast and off the areola. Place your thumb on top of the breast and gently on the tube (don't press hard with your thumb or the tube will point upward). Tickle baby's lower lip with your nipple and wait till the baby opens wide. Position the tube under the baby's upper lip as you pull baby close to the breast so that baby accepts both the breast and the tube. The tip of the baby's nose should just touch the breast. This should help baby accept both your breast and the tube at the same time. If you have difficulty, don't worry as this may take a little practice. Remove baby and start again.

9 It may be necessary to prime the system by squeezing gently on bottle. This will cause the tubing to be filled with supplement and prevent the infant from swallowing air.

10 After use, wind tubing around attachment ring and close system with cover.

CLEANING & STERILIZATION

1. Hospital Sterilization

AUTO-CLAVE

All disassembled parts can be autoclaved to a maximum of 272°F for 3 minutes at 29 PSIG or 250°F for 15 minutes at 15 PSIG.

To avoid damage during sterilization

- Sterilize parts without applying external pressure on parts.

GAS

If gas sterilization is used, residual testing is recommended to determine the proper length of aeration.

2. Home Use

Prior to first use:

1. Wash and rinse parts.
2. Place parts in a 4- to 5-quart saucepan.
3. Fill to within 1 inch of rim with cold water.
4. Cover and boil for 20 minutes.
5. Drain water off, allow parts to cool in open pan.
6. Remove all parts and allow to air-dry.

BOIL

After each use:

1. Fill SNS full of warm soapy water, replace attachment ring and tubing and force water through tubing.
2. Refill with clean water and again force water through tubing.
3. Squeeze the empty bottle a few times to remove water droplets from the tubings.
4. Completely disassemble and clean remaining parts with soapy water.
5. Rinse well and drain dry on a clean towel.
6. Reassemble.
7. It is not necessary to boil parts after each feeding unless your health care professional specifically directs you to do so.
8. If you find mineral deposits building up, try rinsing with vinegar and then thoroughly rinsing out. Alternatively, you can boil with distilled water only.

WASH

RINSE

AIR-DRY

Tube Position

When first using your SNS, you will want to know just how the tubing should be placed. Make sure the baby keeps the tube in the center of its upper lip and not in the corner of its mouth.

Tubing Size

Babies using SNS may need different tubing sizes depending on the exact nature of the problem. The larger the diameter of the tubing, the faster the flow and the less work that the baby has to do. Generally, you will want to begin with the Middle size tubing (see flow regulation). You can judge if you have the correct sized tubing by observing the baby's sucking and swallowing pattern. Some babies may not swallow with every suck or even every other suck. You want to achieve swallowing with every 1-3 suckles. If baby doesn't swallow consistently make sure the tube is positioned properly or try a larger tube. You may have to experiment in order to find the correct size tubing. Your health care professional will probably select the tube size based on the baby's and mother's exact problem. Please read the flow regulation on specific problems for more information. Some general recommendations:

Small Size Suggestions (red valve):

- 1. Baby gaining well but sucking needs to be improved.**
- 2. Thinner fluids such as water. (Water would be used when the baby does not use supplement but needs a better flow to improve sucking.)**

Middle Size Suggestions (white valve):

- 1. Beginning Use**
- 2. Adoptive Nursing**
- 3. Slow Gaining Babies**
- 4. Most Sucking Problems**
- 5. Premature Babies**

Large Size Suggestions (transparent valve):

- 1. Cleft Palate Babies**
- 2. Sucking Problems characterized by a tight mouth.
(Faster flow prevents tightening of mouth.)**
- 3. Severe Hypotonia that results in poor action
by the mouth.**
- 4. Weak, Lazy Sucking Babies**

However, the exact size of tubing needed will depend on several factors, including strength of baby's suck, reason for SNS use, etc. It will probably be necessary for you to experiment with several sizes of tubing. Try the middle size first. If the baby is not swallowing on every third suck, try the larger size. If you aren't sure what swallows sound like, ask your health care professional to help you.

Flow Regulation

The Medela SNS is a very versatile device. It allows the flow to be regulated in a number of different ways. The easiest way to adjust the flow is by the size tubing used (see previous section on tubing size). The larger the tubing the faster the supplement will flow. Three sizes of tubing are included. This will enable you to adjust the flow to the baby's need and the particular problem encountered.

Another way to adjust the flow is by the height of the device. When the bottom of the SNS bottle is above the level of the nipple, the supplement will flow more quickly and will continue to flow into the baby's mouth even when sucking has topped. Such a non-stop flow can be uncomfortable so use extreme caution when the device is above the level of your nipple, especially if the baby seems to loudly gulp or choke. (Most babies will stop and pull away if this should happen.)

Similarly, when the SNS is placed below the level of the nipple, the baby must suck against the flow of gravity. It could be difficult for the baby to suck the supplement in this case. If your health care professional suggests the technique, make sure the baby swallows with every suck or every other suck.

The flow can also be increased by squeezing the bottle in some cases. Place the heel of the hand against bottle while supporting the breast and gently compress the SNS to force more supplement through the device. This could be used when the infant has a very weak suck or needs that extra assistance, such as with the reluctant nurser.

Some mothers have found that warming the supplement in the SNS will increase the flow. It is usually not necessary to do this. However, it is an additional way to increase flow.

The versatility of the Medela SNS allows you to vary the rate of flow: **To summarize**

Increase Flow

1. Use larger tubing
2. Raise device
3. Squeeze bottle
4. Lengthen tube
5. Warm supplement

Decrease flow

1. Smaller size tubing
2. Lower device

SPECIFIC PROBLEMS

Adoptive Mother

It is possible for an adoptive mother to induce lactation. The breasts are quite unique in that if stimulation is applied, a signal is sent to the brain that helps prepare the breasts for nursing. This is true even if the mother is not pregnant or has never been pregnant. It might be easier for the mother who has been pregnant but can also be used for someone who has never been pregnant. The same signal then causes the glands to begin producing milk.

When you first decide that you do want to nurse your adoptive baby you may choose to begin stimulation with the aid of a high quality electric breastpump such as the Medela Lactina® Electric Breastpump. The pump can simulate the baby's sucking stimulation and begin the milk production process. (See appendix for referral for electric pumps.)

If your adoptive plans leave no time for this preparation, don't worry as it is not vital to your success at adoptive nursing. Many successful adoptive breastfeeding mothers have not felt the need to do this, nor have they pumped in between feedings. Consider this only for extra stimulation.

The longer you use your SNS with your adoptive baby, the better the chance that the stimulation will help bring in your milk supply. Most women will not be able to exclusively breastfeed their adoptive infants, so some supplementation is often required especially until 4-6 months or until solid food is introduced and perhaps beyond that. Frequent weight checks with your doctor are very reassuring that baby is doing well. As your milk comes in, your baby will take less from the SNS. Let the baby do this natural weaning from the SNS, and don't arbitrarily cut down on the supplement.

Cleft Palate

Nursing the cleft palate baby can be a most gratifying experience. However, you may encounter some special problems. Because of the infant's difficulty in obtaining a good suction, it may be necessary to use some techniques to increase flow. This includes raising the device, using a larger diameter tubing, and squeezing. See the flow regulation section for more details.

Reluctant Nurser

Some babies have a difficult time learning to breastfeed. In this case, the SNS can be very useful. By attaching the SNS, the infant gets immediate gratification and this may help the baby realize the breast is the sources of milk. Be sure to prime the tube by squeezing the bottle so that the supplement is ready at the end of the tube to reward baby for her first suck. A football hold can offer better control of baby's head if baby is confused and moves head rapidly from side to side. If baby seems frantic, wrap tightly in a receiving blanket to control baby's flailing arms. Return baby to football hold after calming by rocking.

Building Milk Supply

Occasionally mother-baby separation, breastfeeding problems, illness, etc., will cause a temporary reduction in milk supply. The SNS will nourish the baby while at the same time create frequent sucking stimuli which will rapidly increase the milk supply. At first, baby will greedily take the supplement at the breast. At about 48 hours after you begin to use the SNS (depending on the extent of milk supply decrease) you will notice a decrease in the amount of supplement taken. This is an indication that your supply is increasing especially if baby swallows with

every suck. Let baby self wean from the SNS by discontinuing use if baby takes less than 2-3 ounces from the SNS every day. Weight checks with your healthcare provider can be reassuring at this time.

Premature Infant

The SNS can greatly help with breastfeeding as your baby is weaned from gavage feeding in the NICU. Premature babies generally do not have a lot of stamina. If they suck at the breast, premature babies should get some nutrition. The SNS will allow you to help your infant learn to suck at the breast while receiving nutrition that the hospital can monitor. This saves precious energy over breastfeeding alone followed by bottle supplement.

It's best to supplement while nursing, rather than let the infant nurse first and then bring in the SNS. The baby's energy requirements must be taken into account. Short frequent nursing is generally the rule, say at least 10 minutes on one breast, at 1-2 1/2 hour intervals depending on the size of your baby and how well and how long it takes to empty baby's stomach. Premature babies rest and pause a lot during feeding and need a lot of careful monitoring so they actually feed. Pausing is okay and necessary but a gentle nudge may be needed if baby's pausing lasts more than a minute.

Anytime that the infant can suck on a bottle, she can also suck at the breast. With the SNS there is no reason to first feed the infant by bottle and then later transfer baby to the breast. Transferring the baby from bottle to breast runs the risk of nipple confusion/preference.

Slow Gaining Breastfed Baby

If you have been nursing frequently (every 2-3 hours round the clock) and for long feedings (20-60 minutes) and your baby has not gained well, your doctor may suspect illness in the baby or perhaps a sucking problem. The SNS will provide the nutrition needed to your baby so that you can continue to breastfeed, while your healthcare provider investigates the cause of the problem.

If your baby has a sucking in-coordination or sucking problem, the baby will often suck better with the SNS in place for every feeding. Your healthcare provider will probably explain that most sucking problems are temporary and usually resolve when a baby's coordination has a chance to mature. Your healthcare provider may not be able to tell you exactly when your baby will outgrow this in-coordination and when you can discontinue using your SNS. This may take a few weeks to a few months.

Since you will be using your SNS until baby becomes more coordinated, it is advisable to let baby wean herself from the SNS as her suck improves. You should notice an immediate and gratifying weight gain each week. Even though baby is doing well, do not remove the SNS or cut down on supplement at this time. Remember that baby has a sucking problem and can't suck well. Removing the SNS too soon might return the baby to a poor weight gain situation. When baby gains better sucking coordination, baby will, on his/her own, take less supplement. You might try starting the baby at the breast without the SNS. When the suck-swallow rhythm starts slowing down and you no longer hear the swallows, you can then add the SNS to finish the feeding. Remember, do not remove SNS without consulting with your healthcare provider.

Babies with suckling problems may need different tubing sizes depending on the exact nature of the problem

(see flow regulation and tubing size section). Try different sized tubings. You want to be sure that the baby gets some supplement with every 1-3 sucks. If baby doesn't swallow consistently, try a larger tube. Because baby is uncoordinated, she/he may not always swallow with every suck or every other suck.

DECREASING SUPPLEMENTS

One of the goals of using the SNS is to return the mother and baby to as close to a natural nursing situation as possible. The easiest way to do this is to watch baby decrease amount of supplement as your body begins producing more milk and/or baby's suck improves. Several indications are:

Baby gaining weight

— 4-7 ounces per week for breastfed babies
is average

Feelings of milk letdown or leaking

— it can be normal not to notice this and still
produce a sufficient amount of milk

**Baby leaving supplement in SNS but being satisfied
with feeding**

Noticing yourself being able to pump more milk

Feeling more breast fullness

— but not engorgement

You should never put less supplement in SNS than baby can take. Do not attempt to wean baby from SNS by doing this. Let baby wean naturally from the SNS. As you begin to produce more milk, baby will naturally need less and less from the SNS. Let your health care professional determine when you can discontinue using your SNS.

BREASTPUMPS

If you wish to use breastmilk with your SNS, you may want to use an electric breastpump. Most women find this is most convenient and easy because both breasts can be pumped simultaneously. Using a double pump system cuts pumping time in half. In addition, research indicates that double-pumping increases prolactin levels, which helps maintain your milk.

Medela breastpumps can be rented at thousands of independent rental stations throughout the United States and Canada. Call toll free 1-800-TELL YOU (1-800-835-5968) for a rental station in your area.

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©2001 Medela, Inc.
SNS and Supplemental Nursing System are trademarks of Medela, Inc.
Medela is a registered trademark of Medela, Inc.
Medela, Inc.
Medela Canada, Inc.